
ISX-3

产品编号: D52478

CAS: 912789-08-3

分子式: C₁₆H₁₅ClN₄O₂

纯度: 99.32%

InChi Key: InChIKey=DPAAHSSAGAOLBZ-UHFFFAOYSA-N

Smiles: O=C(NCCCN1C=NC=C1)C2=NOC(=C2)C=3C=CC=C(Cl)C3

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 11 mg/mL (33.26 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: ISX-3 is a potent pro-osteogenic and anti-adipogenic agent. ISX-3 can increase the expression of PPAR γ . ISX-3 can be for the research of osteopenia and osteoporosis.