
CVN424

产品编号: D52485

CAS: 2254706-21-1

分子式: C₂₄H₂₉F₂N₅O₃

纯度: 98.8%

InChi Key: InChIKey=HSWVJQBEXRKOBZ-QGZVFWFLSA-N

Smiles: CC(N(CC1)Cc(nc2N[C@H]3COCC3)c1nc2N(CC1)CCC1Oc(ccc(F)c1)c1F)=O

外观: 固体粉末

作用通路: Endocrinology/Hormones;GPCR/G Protein

作用靶点: GPR

溶解性: DMSO :60mg/mL (126.71 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CVN424 is a selective and novel GPR6 Inverse agonist effective in models of Parkinson Disease.