
TAK-041

产品编号: D52508

CAS: 1929519-13-0

分子式: C₁₈H₁₅F₃N₄O₃

纯度: 99.49%

InChi Key: InChIKey=JZGLECLGVQRPPI-NSHDSACASA-N

Smiles: C[C@@H](c(cc1)ccc1OC(F)(F)F)NC(CN1N=Nc(cccc2)c2C1=O)=O

外观: 固体粉末

作用通路: Endocrinology/Hormones; GPCR/G Protein

作用靶点: GPR

溶解性: DMSO: 100 mg/mL (254.89 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TAK-041 is a potent and selective GPR139 agonist.