
Senaparib

产品编号: D52527

CAS: 1401682-78-7

分子式: C₂₄H₂₀F₂N₆O₃

纯度: 99.41%

InChi Key: InChIKey=VBTUJTGLLREMNW-UHFFFAOYSA-N

Smiles: O=C(c(cc(CN(c1cccc(F)c1C(N1)=O)C1=O)cc1)c1F)N(CC1)CCN1c1ncccn1

外观: 固体粉末

作用通路: Chromatin/Epigenetic; DNA Damage/DNA Repair

作用靶点: PARP

溶解性: DMSO:75 mg/ml (156.76 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Senaparib is a novel highly potent and selective oral PARP1/2 inhibitor with strong antitumor activity.