
PSN 375963

产品编号: D52583

CAS: 388575-52-8

分子式: C₁₇H₂₃N₃O

纯度: 99.24%

InChi Key: InChIKey=OAVLEYPTWABFLF-UHFFFAOYSA-N

Smiles: CCCCC(CC1)CCC1c1nc(-c2ccncc2)no1

外观: 固体粉末

作用通路: Endocrinology/Hormones; GPCR/G Protein

作用靶点: GPR

溶解性: DMSO: 4 mg/ml(12.43 mM), Ethanol: 24 mg/ml(74.57 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PSN 375963 hydrochloride is a synthetic agonist of the endogenous ligand for GPR119.