

---

## Fosinopril sodium

产品编号: D52585

CAS: 88889-14-9

分子式: C<sub>30</sub>H<sub>45</sub>NO<sub>7</sub>P·Na

纯度: 99.84%

InChi Key: InChIKey=TVTJZMHAIQQZTL-HREVRLCXS-M

Smiles: [Na+].CCC(=O)OC(OP(=O)(CCCCc1ccccc1)CC(=O)N1C[C@@H](C[C@H]1C([O-])=O)C1CCCCC1)C(C)C

外观: 固体粉末

作用通路: Endocrinology/Hormones

作用靶点: RAAS

溶解性: DMSO:2.5 mM

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Fosinopril sodium is a phosphinic acid-containing angiotensin-converting enzyme inhibitor that is effective in the treatment of hypertension. It is a prodrug that is converted to its active metabolite fosinoprilat.