
Tazarotenic acid

产品编号: D52601

CAS: 118292-41-4

分子式: C₁₉H₁₇NO₂S

纯度: 98.64%

InChi Key: InChIKey=IQIBKLWBVJPOQO-UHFFFAOYSA-N

Smiles: CC1(C)CCSc2ccc(C#Cc3ccc(C(=O)O)cn3)cc21

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 10 mg/mL (30.92 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: AGN-190299 is a retinoid prodrug which is converted to its active form