
PiB

产品编号: D52670

CAS: 566169-93-5

分子式: C₁₄H₁₂N₂O₂S

纯度: 97.02%

InChi Key: InChIKey=ZQAQXZBSGZUUNL-UHFFFAOYSA-N

Smiles: OC1=CC=C2N=C(C3=CC=C(NC)C=C3)SC2=C1

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 60 mg/ml (234.08 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: PiB is a standard for analysis of ¹¹C-PiB. ¹¹C-PiB is a β-amyloid PET imaging tracer for Alzheimer's disease diagnosis.