
MP7

产品编号: D52672

CAS: 1001409-50-2

分子式: C₂₈H₂₂F₂N₄O₄

纯度: 99.66%

InChi Key: InChIKey=GCWCGSPBENFEPE-VWLOTQADSA-N

Smiles: O=C(C1=CC=CN(Cc(cc2)cc(F)c2F)C1=O)N[C@@H](COc(cc1)cc(N2)c1NC2=O)c1cccc1

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PDK

溶解性: DMSO: 100 mg/mL (193.61 mM); H₂O: < 0.1 mg/mL (insoluble)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: MP7 is a inhibitor of phosphoinositide-dependent kinase-1 (PDK1).