
M1001

产品编号: D52676

CAS: 874590-32-6

分子式: C₁₇H₁₇N₃O₂S

纯度: 100%

InChi Key: InChIKey=LONXCUOJVJLTIP-UHFFFAOYSA-N

Smiles: O=S1(=O)N=C(NC2=CC=CC=C2N2CCCC2)C2=CC=CC=C12

外观: 固体粉末

作用通路: Angiogenesis; Chromatin/Epigenetic; Metabolism

作用靶点: HIF; HIF/HIF Prolyl-Hydroxylase

溶解性: DMSO:250 mg/mL (763.59 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: M1001 is a HIF-2 α agonist.