

---

## TSHR antagonist S37

产品编号: D52685

CAS: 1217616-61-9

分子式: C<sub>25</sub>H<sub>20</sub>N<sub>2</sub>O<sub>3</sub>S<sub>2</sub>

纯度: 98.35%

InChi Key: InChIKey=YGFJFPYQZCZNIH-UHFFFAOYSA-N

Smiles: O=C(N1)SC(C2C3=CC=CC=C3)=C1SC(C2C4C56)C(C4)C6C(N(C7=CC=CC=C7)C5=O)=O

外观: 固体粉末

作用通路: GPCR/G Protein

作用靶点: TSH Receptor

溶解性: DMSO: 65 mg/ml (141.13 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TSHR antagonist S37 is a selective and competitive antagonist of the thyrotropin receptor (TSHR).