
trans-isrib A17

产品编号: D52689

CAS: 1628478-15-8

分子式: C₂₂H₂₂Cl₂F₂N₂O₄

纯度: 98.17%

InChi Key: InChIKey=BLLMXZJTGLCEBP-HDJSIYSDSA-N

Smiles: O(CC(N[C@H]1CC[C@H](NC(COC2=CC(F)=C(Cl)C=C2)=O)CC1)=O)C3=CC(F)=C(Cl)C=C3

外观: 固体粉末

作用通路: Others

作用靶点: Others

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: trans-isrib A17 is an potent inhibitor of the integrated stress response with EC₅₀ of 0.6 nM.