
Adoprazine

产品编号: D52716

CAS: 222551-17-9

分子式: C₂₄H₂₄N₃O₂

纯度: 99.44%

InChi Key: InChIKey=IUVSEUFHPNITEQ-UHFFFAOYSA-N

Smiles: FC1=CC=C(C=C1)C=2C=NC=C(C2)CN3CCN(C=4C=CC=C5OCCOC54)CC3

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 45mg/mL (111mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Adoprazine is a potential atypical antipsychotic bearing potent D2 receptor antagonist and 5-HT_{1A} receptor agonist properties.