
Intepirdine

产品编号: D52800

CAS: 607742-69-8

分子式: C₁₉H₁₉N₃O₂S

纯度: 95.19%

InChi Key: InChIKey=JJZFWROHYSMCMU-UHFFFAOYSA-N

Smiles: c1(cccc2c1ncc(c2)S(=O)(=O)c1ccccc1)N1CCNCC1

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: 5-HT Receptor

溶解性: DMSO:72 mg/mL (203.7 mM);Ethanol:2 mg/mL (5.65 mM);Water:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: SB742457 is a highly selective 5-HT₆ receptor antagonist, used in trials studying the treatment of Alzheimer's Disease.