
IPI-3063

产品编号: D52805

CAS: 1425043-73-7

分子式: C₂₅H₂₅N₇O₂

纯度: 99.77%

InChi Key: InChIKey=OBHAYOJCPNWKBL-HNNXBMFYSA-N

Smiles: C[C@@H](C1=CC2=C(C(=CC=C2)C3=CN(C(=O)C=C3)C(=O)N1C(C)C)NC4=NC=NC(=C4C#N)N

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling

作用靶点: PI3K

溶解性: DMSO:83.33 mg/mL (182.94 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: IPI-3063 is an effective and selective inhibitor of PI3K p110 δ (IC₅₀: 2.5 $\pm$ 1.2 nM).