
TC-DAPK 6

产品编号: D52823

CAS: 315694-89-4

分子式: C₁₇H₁₂N₂O₂

纯度: 98.13%

InChi Key: InChIKey=GFGMISOSPOPSHN-UHFFFAOYSA-N

Smiles: O=C1OC(=NC1=CC=2C=NC=CC2)C=CC=3C=CC=CC3

外观: 固体粉末

作用通路: Apoptosis

作用靶点: DAPK

溶解性: DMSO:16.8 mg/mL (60.8 mM; Need ultrasonic and warming)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TC-DAPK 6 is an ATP-competitive inhibitor of Death-associated protein kinase (DAPK) with IC₅₀s of 69 and 225 nM for DAPK1 and DAPK3.