
WH-4-025

产品编号: D52833

CAS: 1876463-35-2

分子式: C₃₉H₃₈F₃N₇O₅

纯度: 98.7%

InChi Key: InChIKey=MWWDOHSKUDECfq-UHFFFAOYSA-N

Smiles: O=C(OC1=CC(C(NC2=CC=CC(F)(F)F)=C2)=O)=CC=C1C)N(C3=CC=C(OC)C=C3OC)C4=NC(NC5=CC=C(N6CCN(C)CC6)C=C5)=NC=C4

外观: 固体粉末

作用通路: PI3K/Akt/mTOR signaling; Chromatin/Epigenetic

作用靶点: AMPK; SIK

溶解性: DMSO:125 mg/mL (164.47 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: WH-4-025 is a Salt-inducible kinase (SIK) inhibitor.