
BE1218

产品编号: D52834

CAS: 2893967-40-1

分子式: C₃₀H₃₀FNO₄S₂

纯度: 99.7%

InChi Key: InChIKey=YZKNHPOMZOAMSW-UHFFFAOYSA-N

Smiles: O=S(=O)(C=1C(=CC(=CC1C)C)C)N(CC2=CC=C(C=C2)C=3C=CC=C(C3)S(=O)(=O)C)CC=4C=CC=CC4F

外观: 固体粉末

作用通路: Metabolism

作用靶点: Liver X Receptor

溶解性: DMSO: 5.52 mg/mL (10 mM), Need ultrasonic

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: BE1218 is a liver X receptor (LXR) inverse agonist active on LXR α and LXR β with an IC₅₀ of 9 nM and 7 nM, respectively.