
KGP94

产品编号: D52840

CAS: 1131456-28-4

分子式: C₁₄H₁₂BrN₃OS

纯度: 98.52%

InChi Key: InChIKey=ZDBKSZKTCPOBFR-UHFFFAOYSA-N

Smiles: NC(NN=C(c1cc(O)ccc1)c1cccc(Br)c1)=S

外观: 固体粉末

作用通路: Proteases/Proteasome

作用靶点: Cysteine Protease

溶解性: DMSO: 60 mg/ml (171.32 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: KGP94 is a potent, selective, and competitive inhibitor of the lysosomal endopeptidase enzyme.