
AZD4320

产品编号: D52854

CAS: 1357576-48-7

分子式: C₄₅H₄₈ClF₃N₄O₇S₃

纯度: 99.71%

InChi Key: InChIKey=UNEJSHNDABUZYUJNHCCGESAN

Smiles: CN(CCO)CC[C@H](CSc1ccccc1)Nc1ccc(cc1S(=O)(=O)C(F)(F)F)S(=O)(=O)NC(=O)c1ccc(cc1)N1C
CC(CC1)[C@@H](O)c1ccccc1-c1ccc(Cl)cc1

外观: 固体粉末

作用通路: Apoptosis

作用靶点: BCL

溶解性: DMSO:125 mg/mL (132.20 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: AZD4320 is an inhibitor of BH3-mimicking dual BCL2/BCLxL. For KPUM-MS3, KPUM-UH1, and STR-428 cells, the IC₅₀s values are 26 nM, 17 nM, and 170 nM , respectively.