
AS601245.2TFA (345987-15-7 free base)

产品编号: D52857

CAS: 345987-16-8

分子式: C₂₄H₁₈F₆N₆O₄S

纯度: 98%

InChi Key: InChIKey=CIMRIDPPEONTRF-UHFFFAOYSA-N

Smiles: OC(=O)C(F)(F)F.OC(=O)C(F)(F)F.N#CC(c1nc2ccccc2s1)c1ccnc(NCCc2cccnc2)n1

外观: 固体粉末

作用通路: MAPK

作用靶点: JNK

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: AS601245.2TFA is a cell-permeable Inhibitor of JNK (IC₅₀s of 150, 220, and 70 nM for hJNK1, hJNK2, and hJNK3, respectively).