
SPD304 dihydrochloride

产品编号: D52870

CAS: 1049741-03-8

分子式: C₃₂H₃₄Cl₂F₃N₃O₂

纯度: 98%

InChi Key: InChIKey=GOZMBJCYMQQACI-UHFFFAOYSA-N

Smiles: O=C1C2=C(OC=C1CN(CCN(CC3=CN(C4=C3C=CC=C4)C5=CC=CC(C(F)(F)F)=C5)C)C=C(C)C(C)=C2.Cl.Cl

外观: 固体粉末

作用通路: Apoptosis

作用靶点: TNF

溶解性: DMSO: 6.2 mg/mL (10 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: SPD304 dihydrochloride is a selective inhibitor of TNF- α with an IC₅₀ of 22 μ M. SPD304 dihydrochloride promotes dissociation of TNF trimers.