
TAS-103 dihydrochloride

产品编号: D52888

CAS: 174634-09-4

分子式: C₂₀H₂₁Cl₂N₃O₂

纯度: 98%

InChi Key: InChIKey=HAYAYGFVSIWSGQ-UHFFFAOYSA-N

Smiles: Cl.Cl.CN(C)CCNc1nc2cc(O)ccc2c2-c3ccccc3C(=O)c12

外观: 固体粉末

作用通路: DNA Damage/DNA Repair

作用靶点: Topoisomerase

溶解性: H₂O:50 mg/mL (123.06 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TAS-103 (dihydrochloride) is a novel anticancer agent targeting both topoisomerase (Topo) I and Topo II.