
CASP3 Activator 1541

产品编号: D52892

CAS: 1100353-03-4

分子式: C₂₄H₁₇N₃O₄

纯度: 99.98%

InChi Key: InChiKey=PQLVYQVWNPGJML-UHFFFAOYSA-N

Smiles: COc1cccc2cc(C(=O)Nc3cccc(c3)-c3cn4cccc4n3)c(=O)oc12

外观: 固体粉末

作用通路: Apoptosis; Proteases/Proteasome

作用靶点: Caspase

溶解性: DMSO: 4.11 mg/ml (10 mM, Need ultrasonic); Water: Insoluble

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CASP3 Activator 1541 is an activator of the proenzyme forms of caspase-3, and induces cell death.