
Bavisant

产品编号: D52920

CAS: 929622-08-2

分子式: C₁₉H₂₇N₃O₂

纯度: 96.06%

InChi Key: InChIKey=BGBVSGSIXIIREO-UHFFFAOYSA-N

Smiles: O=C(c1ccc(CN2CCOCC2)cc1)N(CC1)CCN1C1CC1

外观: 固体粉末

作用通路: GPCR/G Protein;Immunology/Inflammation;Neuroscience

作用靶点: Histamine Receptor

溶解性: DMSO: 13 mg/mL (39.46 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Bavisant is a selective and orally active Human H₃ receptor antagonist. Bavisant can be used in research on mechanisms of action, involving wakefulness and cognition and the treatment for ADHD.