
p38 MAPK Inhibitor

产品编号: D52940

CAS: 219138-24-6

分子式: C₂₀H₁₃ClFN₃O

纯度: 99.24%

InChi Key: InChIKey=DZFBYHUKZSRPHU-UHFFFAOYSA-N

Smiles: O=C1C(C=2C=CC(F)=CC2)=C(NN1C3=CC=C(Cl)C=C3)C=4C=CN=CC4

外观: 固体粉末

作用通路: MAPK

作用靶点: p38 MAPK

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: p38 MAPK inhibitor is a potent inhibitor of p38 MAP kinase (IC₅₀ = 35 nM). It inhibits senescence induced by the oncogene RAS.