
VUF 5574

产品编号: D52950

CAS: 280570-45-8

分子式: C₂₁H₁₇N₅O₂

纯度: 97.63%

InChi Key: InChIKey=YRAFEJSZTVWUMD-UHFFFAOYSA-N

Smiles: O=C(NC1=NC(=NC=CC=CC21)C3=CN=CC=C3)NC=CC=CC4OC

外观: 固体粉末

作用通路: GPCR/G Protein;Neuroscience

作用靶点: Adenosine Receptor

溶解性: DMSO: Insoluble

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: VUF 5574 is a selective antagonist of adenosine A₃ receptor with a K_i of 4.03 nM for the recombinant human receptor.