
CNT2 inhibitor-1

产品编号: D52972

CAS: 880155-70-4

分子式: C₂₃H₂₄N₆O₄

纯度: 99.37%

InChi Key: InChiKey=YJLIKUSWRSEPSM-WGQQHEPDSA-N

Smiles: Nc1ncnc2n([C@@H]3O[C@H](CO)[C@@H](O)[C@H]3O)c(NCc3ccc(cc3)-c3ccccc3)nc12

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 70 mg/mL (156.1 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: CNT2 inhibitor-1 is a potent inhibitor of concentrative nucleoside transporter 2 (CNT2; IC₅₀: 640 nM for hCNT2).