
2-Iodohippuric Acid

产品编号: D52986

CAS: 147-58-0

分子式: C₉H₈INO₃

纯度: 95%

InChi Key: InChIKey=CORFWQGVBFZHF-UHFFFAOYSA-N

Smiles: OC(=O)CNC(=O)c1ccccc1I

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 55 mg/mL (180.29 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: 2-Iodohippuric Acid is an organic compound that can be used to measure renal clearance for renal function studies, and for the diagnosis of certain diseases, such as liver and kidney disease. In addition, it is used in the synthesis of certain drugs.