

---

## WNK463

产品编号: D53002

CAS: 2012607-27-9

分子式: C<sub>21</sub>H<sub>24</sub>F<sub>3</sub>N<sub>7</sub>O<sub>2</sub>

纯度: 100%

InChi Key: InChIKey=HWSHOMMVLGBIDN-UHFFFAOYSA-N

Smiles: O=C(c1cncn1C1CCN(c2ncc(c3nnc(C(F)(F)F)O3)cc2)CC1)NC(C)(C)C

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint; Metabolism

作用靶点: Serine/threonin kinase

溶解性: DMSO:85 mg/mL (183.4 mM);Ethanol:85 mg/mL (183.4 mM);Water:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: WNK463 is a pan-WNK-kinase inhibitor. It effectively inhibits the in vitro kinase activity of all four WNK family members (WNK1/2/3/4).