
TCPOBOP

产品编号: D53043

CAS: 76150-91-9

分子式: C₁₆H₈Cl₄N₂O₂

纯度: 98.08%

InChi Key: InChIKey=BAFKRPOFIYPKBQ-UHFFFAOYSA-N

Smiles: ClC1=CN=C(OC2=CC=C(OC3=NC=C(Cl)C=C3Cl)C=C2)C(Cl)=C1

外观: 固体粉末

作用通路: Apoptosis

作用靶点: BCL

溶解性: DMSO:25 mg/mL (62.18 mM; ultrasonic and warming and heat to 50°C)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TCPOBOP is a constitutive androstane receptor agonist that induces robust hepatocyte proliferation and hepatomegaly. TCPOBOP decreases Fas-induced murine liver injury by altering Bcl-2 proteins.