

---

## RORyt Inverse agonist 6

产品编号: D53058

CAS: 1887161-80-9

分子式: C<sub>28</sub>H<sub>29</sub>ClN<sub>6</sub>O<sub>5</sub>

纯度: 98%

InChi Key: InChIKey=QUTSLTZISBFEOJ-XMMPXPASA-N

Smiles: C(C)(C)N1C=2C(C(=O)N(CC3CC3)C1=O)=CC(NC(=O)N4C[C@H](C(NC5=CC(Cl)=C(C#N)C=C5)=O)OCC4)=CC2

外观: 固体粉末

作用通路: Metabolism

作用靶点: ROR

溶解性: DMSO : 45 mg/mL (79.64 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: RORyt Inverse agonist 6 is an agonist of RORyt inverse. RORyt Inverse agonist 6 can be used in research on Th17-driven autoimmune diseases.