
MSC-4106

产品编号: D53098

CAS: 2738542-58-8

分子式: C₁₈H₁₂F₃N₃O₂

纯度: 100%

InChi Key: InChIKey=HUSMWXZHQVTDRU-UHFFFAOYSA-N

Smiles: O=C(O)C=1C=CC2=C(C1)C3=CN(N=C3N2C4=CC=C(C=C4)C(F)(F)F)C

外观: 固体粉末

作用通路: Stem Cells

作用靶点: YAP

溶解性: DMSO: 225mg/mL (626.2mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: MSC-4106 is an orally active and potent inhibitor of YAP/TAZ-TEAD. MSC-4106 inhibits TEAD1 or TEAD3 auto-palmitoylation and shows inhibitory effect on NCI-H226 tumor xenograft model.