

---

## DN401

产品编号: D53100

CAS: 2135749-60-7

分子式: C<sub>18</sub>H<sub>19</sub>N<sub>3</sub>O<sub>2</sub>

纯度: 100%

InChi Key: InChIKey=WWQDEXGFVSTCX-UHFFFAOYSA-N

Smiles: N=1N=C(C)CC=2C=C(OC)C(OC)=CC2C1C=3C=CC(N)=CC3

外观: 固体粉末

作用通路: Cytoskeletal Signaling;Metabolism

作用靶点: HSP

溶解性: DMSO: 3.83 mg/mL (10 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: DN401 is a potent TRAP1 and Hsp90 Inhibitor( IC<sub>50</sub> = 79 nM for TRAP1, IC<sub>50</sub> = 698 nM for Hsp90) with potent anticancer activity.