
Ningetinib Tosylate

产品编号: D53106

CAS: 1394820-77-9

分子式: C₃₈H₃₇FN₄O₈S

纯度: 99.77%

InChi Key: InChIKey=NAUYISNCVNUUDK-UHFFFAOYSA-N

Smiles: O=C(C1=C(C)N(C)N(C2=CC=CC=C2)C1=O)NC3=CC=C(OC4=CC=NC5=CC(OCC(C)(O)C)=CC=C45)C(F)=C3.CC6=CC=C(S(=O)(O)=O)C=C6

外观: 固体粉末

作用通路: Angiogenesis; Tyrosine Kinase/Adaptors

作用靶点: c-Met/HGFR; TAM Receptor; VEGFR

溶解性: DMSO:8 mg/mL (10.98 mM; Need ultrasonic and warming)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Ningetinib Tosylate is an orally bioavailable tyrosine kinase inhibitor with IC₅₀s of <1.0, 1.9 and 6.7 nM for Axl, VEGFR2, and c-Met, respectively.