
zapnometinib

产品编号: D53121

CAS: 303175-44-2

分子式: C₁₃H₇ClF₂INO₂

纯度: 99.66%

InChi Key: InChIKey=XCNBGWKQXRQKSA-UHFFFAOYSA-N

Smiles: OC(c(c(Nc(ccc(I)c1)c1Cl)c1F)ccc1F)=O

外观: 固体粉末

作用通路: MAPK; Microbiology/virology

作用靶点: Antibacterial; Influenza Virus; MEK

溶解性: DMSO : 60 mg/mL (146.50 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: zapnometinib is a MEK inhibitor.