
Befotertinib

产品编号: D53135

CAS: 1835667-63-4

分子式: C₂₉H₃₂F₃N₇O₂

纯度: 99.4%

InChi Key: InChIKey=USOCZVZOXKTJTI-UHFFFAOYSA-N

Smiles: C=CC(=O)NC1=CC(NC2=NC=CC(C3=CN(CC(F)(F)F)C4=C3C=CC=C4)=N2)=C(OC)C=C1N(C)CCN(C)C

外观: 固体粉末

作用通路: Angiogenesis; JAK/STAT signaling; Tyrosine Kinase/Adaptors

作用靶点: EGFR

溶解性: DMSO : 60 mg/mL (105.71 mM; ultrasonic and warming and heat to 60°C)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Befotertinib is an inhibitor of EGFR tyrosine kinase and can be used for studies about EGFR T790M-positive non-small cell lung cancer.