
Phosphoglucomutase from rabbit muscle

产品编号: D10238

CAS: 9001-81-4

保存条件: Store at 2-8°C.

产品介绍: CAS: 9001-81-4

MDL: MFCD00131886

grade: BR, 100u/mg protein

form: Ammonium sulfate suspension, Crystalline suspension in 3.2 M (NH₄)₂SO₄, pH 6.0
containing 0.01% EDTA

foreign activity:

lactic dehydrogenase ≤0.5%

phosphoglucose isomerase ≤0.01%

pyruvate kinase ≤0.05%

shipped in: wet ice

Unit Definition: One unit will convert 1.0 μmole of α-D-glucose 1-phosphate to α-D-glucose 6-phosphate per min at pH 7.4 at 30°C.