

---

## ElG121 Rabbit pAb

Catalog Number: bs-6533R

Target Protein: ElG121

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500),  
ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Pig, Cow, Dog)

Predicted MW: 107 kDa

Subcellular Cell membrane

Locations:

Entrez Gene: 57535

Swiss Prot: Q6UXG2

Source: KLH conjugated synthetic peptide derived from human Estrogen-induced gene 121 protein:  
841-940/1013.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Estrogen-induced gene 121 protein may play a role as a marker of hyperestrogenic state and estrogen-related type I endometrial carcinoma. It belongs to the UPF0577 family.