
Z-HA155

产品编号: D53780

CAS: 1312201-00-5

分子式: C₂₄H₁₉BFNO₅S

纯度: 98.28%

InChi Key: InChIKey=BRWUZCBSWABPMR-XKZIYDEJSA-N

Smiles: C(=C\1/C(=O)N(CC2=CC=C(F)C=C2)C(=O)S1)\C3=CC=C(OCC4=CC=C(B(O)O)C=C4)C=C3

外观: 固体粉末

作用通路: Metabolism

作用靶点: PDE

溶解性: DMSO: 8 mg/mL (17.27 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: HA155 selectively and potently inhibits autotaxin with an IC₅₀ of 5.7 nM.