
Tirofiban

产品编号: D53804

CAS: 144494-65-5

分子式: C₂₂H₃₆N₂O₅S

纯度: 99.87%

InChi Key: InChIKey=COKMIXFXJJXBQG-NRFXNRHFSA-N

Smiles: CCCCS(=O)(=O)N[C@@H](CC1=CC=C(C=C1)OCCCCC2CCNCC2)C(=O)O

外观: 固体粉末

作用通路: Cytoskeletal Signaling

作用靶点: Integrin

溶解性: DMSO:<1 mg/mL; Ethanol:<1 mg/mL; Water:<1 mg/mL

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Tirofiban (MK-383) is a selective platelet GPIIb/IIIa antagonist which inhibits platelet aggregation with IC₅₀ of 9 nM.