
Talinolol

产品编号: D53805

CAS: 57460-41-0

分子式: C₂₀H₃₃N₃O₃

纯度: 99.15%

InChi Key: InChIKey=MXFWWQICDIZSOA-UHFFFAOYSA-N

Smiles: CC(C)(C)NCC(COc(cc1)ccc1NC(NC1CCCCC1)=O)O

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: Adrenergic Receptor

溶解性: DMSO: 14 mg/mL (38.51 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: (±)-Talinolol is a selective β₁ adrenoceptor antagonist well known for its cardioprotective and antihypertensive activity.