
TCS 2210

产品编号: D53806

CAS: 1201916-31-5

分子式: C₁₈H₁₇N₃O₃

纯度: 97.42%

InChi Key: InChIKey=ZOFZKGFBKISZRI-UHFFFAOYSA-N

Smiles: O=C(NO)C=1C=CC=2NC(=O)C(=NC2C1)CCCC=3C=CC=CC3

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: 4.8 mg/mL (14.84 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TCS 2210 is an inducer of neuronal differentiation in mesenchymal stem cells and increases the neuronal markers β -III tubulin expression and neuron-specific enolase expression.