
TOPS

产品编号: D53813

CAS: 40567-80-4

分子式: C₁₂H₁₈NNaO₃S

纯度: 99.71%

InChi Key: InChIKey=NJZLCEFCAHNYIR-UHFFFAOYSA-M

Smiles: O=S(CCCN(CC)C1=CC=CC(C)=C1)([O-])=O.[Na+]

外观: 固体粉末

作用通路: Others

作用靶点: Others

溶解性: DMSO: ≥ 42 mg/mL (150.36 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: TOPS is a Trinder's reagent and is widely used in diagnostic tests and biochemical tests.