
MF-094

产品编号: D53827

CAS: 2241025-68-1

分子式: C₃₀H₃₇N₃O₄S

纯度: 98.46%

InChi Key: InChiKey=DQXORJHBZHLLTE-SANMLTNESA-N

Smiles: CC(C)(C)NS(c1ccc2cccc1c2NC([C@H](Cc1ccccc1)NC(C1CCCCC1)=O)=O)(=O)=O

外观: 固体粉末

作用通路: Ubiquitination

作用靶点: DUB

溶解性: DMSO: 115 mg/mL (214.67 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: MF-094 is a potent and selective inhibitor of USP30 with IC₅₀ of 120 nM.