
EML 425

产品编号: D53876

CAS: 1675821-32-5

分子式: C₂₇H₂₄N₂O₄

纯度: 100%

InChi Key: InChIKey=LUGQBJDYUPNAQQ-UHFFFAOYSA-N

Smiles: CC1=CC(O)=CC(C)=C1C=C1C(=O)N(CC2=CC=CC=C2)C(=O)N(CC2=CC=CC=C2)C1=O

外观: 固体粉末

作用通路: Chromatin/Epigenetic

作用靶点: Epigenetic Reader Domain; Histone Acetyltransferase

溶解性: DMSO:250 mg/mL (567.55 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: EML 425 is a potent and selective inhibitor of CREB binding protein (CBP)/p300 (IC₅₀s: 2.9 and 1.1 μM, respectively).