

---

## Pizotifen

产品编号: D53877

CAS: 15574-96-6

分子式: C<sub>19</sub>H<sub>21</sub>N

纯度: 99.91%

InChi Key: InChiKey=FIADGNVRKBPQEU-UHFFFAOYSA-N

Smiles: CN1CC(C(CC1)=C1)c2ccsc2CCc2ccccc12

外观: 固体粉末

作用通路: GPCR/G Protein; Neuroscience

作用靶点: 5-HT Receptor

溶解性: H<sub>2</sub>O: < 0.1 mg/mL (insoluble); DMSO: 20 mg/mL (67.70 mM; Need ultrasonic)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Pizotifen is a strong serotonin and tryptamine antagonist, with weak anticholinergic, anti-histaminic and anti-kinin effects. It also possesses sedative and appetite-stimulating properties. Pizotifen hydrochloride is an active ingredient in Sandomigran, which is used for the prophylactic management of migraine. Sandomigran is available in number of countries excluding the United States.