
Dapansutrile

产品编号: D53897

CAS: 54863-37-5

分子式: C₄H₇NO₂S

纯度: 98%

InChi Key: InChIKey=LQFRYKBDZNPJSW-UHFFFAOYSA-N

Smiles: CS(=O)(=O)CCC#N

外观: 固体粉末

作用通路: Immunology/Inflammation; NF-κB

作用靶点: NOD; NOD-like Receptor (NLR)

溶解性: DMSO:120mg/ml (901.1 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Dapansutrile is a potent, selective and orally active NLRP3 inflammasome inhibitor, with Anti-inflammatory and analgesic activity.