
COH34 analog 1

产品编号: D53909

CAS: 314248-24-3

分子式: C₁₈H₁₅NO₃S

纯度: 98%

InChi Key: InChIKey=WCYGPLQLYOSVIH-UHFFFAOYSA-N

Smiles: O=S(=O)(N=CC1=C(O)C=CC=2C=CC=CC21)C3=CC=C(C=C3)C

外观: 固体粉末

作用通路: Others

作用靶点: Others

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: COH34 analog 1 is a COH34 oxidation analogue that can be used to synthesize active compounds.