
Upamostat

产品编号: D53918

CAS: 590368-25-5

分子式: C₃₂H₄₇N₅O₆S

纯度: 98.8%

InChi Key: InChIKey=HUASEDVYRABWCV-NDEPHWFRSA-N

Smiles: CCOC(=O)N1CCN(CC1)C(=O)[C@H](CC2=CC(=CC=C2)/C(=N/O)/N)NS(=O)(=O)C3=C(C=C(C=C3C(C)C)C(C)C)C(C)C

外观: 固体粉末

作用通路: Cell Cycle/Checkpoint; Metabolism; Proteases/Proteasome

作用靶点: PAI-1; Serine Protease; Serine/threonin kinase

溶解性: DMSO:250 mg/mL (396.95 mM)

保存条件: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

产品介绍: Upamostat is an inhibitor of a serine protease and is a urokinase plasminogen activator (uPA) inhibitor.